Digitales

Fieberthermometer

DE

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Art. Nr. 9810102

Hersteller:

Kjümp Health Co., Ltd.

No. 55 Wu Kong Str. Rd. New Taipei Industrial Park, New Taipei City 24890, Taiwan

https://www.polygreen.com

EC REP:

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info@bresser.de

Please read carefully before use

DEVICE DESCRIPTION AND INTENDED USE

This digital thermometer is intended to measure human body temperature for children and adults; it has a display and can be used in household environment.

CONTRAINDICATIONS

- Do not place the thermometer on scarred tissue or tissue compromised by skin disorders.
- Do not place the thermometer probe on skin exposed to direct sunlight, fireplace heat, cold compress therapies, air conditioner flow, etc.
- Do not use **orally**: 1. Infant and children. 2. Unconscious or mouth-breathing patients. 3. Inflammation or surgery of mouth. 4. Persistent frequent coughing. 5. Weak or old patients who can't close the mouth well. 6. Patient with seizure disorder. 7. After drinking hot or cold liquids.
- Do not use **axillary**: 1. Skin disease.
- Do not use **rectal**: 1. Patients with surgical operation in the rectum or perineal region. 2. Disease or inflammation of rectum. 3. Diarrhea or hemorrhoid.

FEATURES

- Easy to read digital LCD. 2. Quick measurement, the minimum measuring time of oral and rectal is around 10 seconds, axillary is around 30 seconds with proper use. 3. Memory display for last temperature taken. 4. Compact, accurate and reasonable device. 5. Very sensitive to body temperature. 6. If the thermometer is inadvertently left with the battery after temperature stabilization, it will automatically shut off in 3 minutes. 7. Small lightweight unit. 8. When low battery indication shows, replace the old battery with a new battery. 9. The unit is water resistant.

PART DESCRIPTION

SPECIFICATIONS

Measuring Range	32.0~43.9 °C
Resolution	0.1 °C
Accuracy	35.5 ~ -42.0 °C ±0.1 °C; Other ±0.2 °C
Memory	Last memory
Battery (included)	1x LR41 Battery
Battery Life	More than 5000 measurements
Operating Condition	16.0 ~ -40.0 °C, R.H. of 15~95% non-condensing
Storage Condition	-20.0 ~ -50.0 °C, R.H. of 15~95% non-condensing
Dimension	150.2mm x 23.2mm x 19.0mm
Weight	Approx. 15g included battery
Certification	ISO 13485
Product Classification	Type BF equipment A

EXPLANATION OF SYMBOLS

	Memory Indication
	Follow operating instructions
	Manufacturer
	Authorized Representative in the European Community.
	Warning and Caution
	MMYY: Lot Number composed by month(M) and year(Y) of manufacture
	The CE marking with the Registration Number of the Notified Body. This denotes the compliance of European Medical Device Directive 93/42/EEC.
	Type BF: Protection against electric shock, type of applied part, body floating.

2. Protected against solid foreign objects of 12.5mm and greater.
2. Dripping water when tilted at an angle.
- Enclosures classified according to degree of protection against ingress of water and particulate matter (IPN12) as equipment per IEC 60529 according to classification in 6.3 (Table D.3, code 2).
- WEEE (Waste Electrical and Electronic Equipment Directive, 2012/19/EU): The device could be send back to the manufacturer for recycling or proper disposal after the expected life time. Alternatively the device shall be disposed in accordance with the national laws.
- Stop using the thermometer if it operates erratically or if the display malfunctions.
- Keep out of reach of unsupervised children. Do not let children play with the unit. TEMPERATURE SHOULD BE TAKEN UNDER ADULT SUPERVISION.
- Open the thermometer before using. Contact a physician immediately should this occur.
- Do not attempt to disassemble the device except to replace the battery.
- Do not use if the display malfunctions or it operates irregularly.
- If the thermometer has been stored at below-freezing temperatures, allow the thermometer to warm naturally to room temperature before using.
- Do not put the device in dishwasher.
- Performance may be degraded.
- RoHS symbol: This denotes the compliance of Restriction of Hazardous Substances. Directive.
- RoHS symbol: This denotes the compliance of Restriction of Hazardous Substances. Directive.
- REACH symbol: Registration, Evaluation, Authorization and Restriction of Chemicals (REACH) is European Union regulation.
- Measuring site is oral, axillary and rectal.
- Medical device

PRECAUTION

- Do not pull the probe. Instead, clean the unit by wiping it with a dry cloth and disinfect the probe with isopropyl (rubbing) alcohol. See cleaning instructions for details.
- Do not drop the thermometer or expose it to heavy shock.
- Do not bite or swallow the probe and the battery cover.
- Do not let thermometer come into contact with any flammable or chemical solvents
- After using, clean the thermometer and place it into the carrier case.
- Do not store the unit under direct sunlight, at a high temperature, in high humidity or dirt. Performance may be degraded.
- Stop using the thermometer if it operates erratically or if the display malfunctions.
- Keep out of reach of unsupervised children. Do not let children play with the unit. TEMPERATURE SHOULD BE TAKEN UNDER ADULT SUPERVISION.
- Open the thermometer before using. Contact a physician immediately should this occur.
- Do not attempt to disassemble the device except to replace the battery.
- Do not use if the display malfunctions or it operates irregularly.
- If the thermometer has been stored at below-freezing temperatures, allow the thermometer to warm naturally to room temperature before using.
- Do not put the device in dishwasher.
- Performance may be degraded.
- RoHS symbol: This denotes the compliance of Restriction of Hazardous Substances. Directive.
- REACH symbol: Registration, Evaluation, Authorization and Restriction of Chemicals (REACH) is European Union regulation.
- Measuring site is oral, axillary and rectal.
- Medical device

- For battery disposal, please consider the following:
- Thermometers contain the battery which can do some measurements and may not offer enough measurements as a new battery.
- If you don't use the thermometer regularly, be sure to remove the battery after each use to avoid the battery leakage or potential damage.

HOW TO USE

- Disinfect the probe with isopropyl (rubbing) alcohol before using.
- Press and release the ON/OFF button, the display read  confirming that the unit is functioning properly.
- Then the display will show L °C with "C" flashing.
- Place the probe in the measurement position.
- Once the degree sign "C" stops flashing, the measurement is finished and the temperature is displayed. An alarm signal will also sound for approximately 8 seconds after "C" stops flashing.
- The unit will automatically turn off in 3 minutes (approx.). However, to prolong battery life, it is best to press the ON/OFF button to turn off once the temperature has been indicated.

Special Features

- Last-Temperature Memory Recall:** When the unit is off, press and hold the ON/OFF button around 3 seconds, the last memory will be shown. Release the ON/OFF button, the unit will be ready for measuring.
- For Models with Fever Alarm:** If the temperature registers higher than 37.5 °C, the fever alarm (a quick beeping sound which occurs continuously for about 8 seconds) will sound when the peak temperature has registered. Also, the unit will keep flashing until the thermometer is turned off.
- NOTE:** Always disinfect the thermometer before using.
- NOTE:** Drinking or eating prior to use may affect reading; oral measurement should be taken 30 minutes after drinking or eating.
- NOTE:** In step 4, you may put a sanitary plastic probe cover on the sensor tip to help to prevent cross-infection and contamination. Use of the probe cover may result in a 0.1 °C difference from actual temperature. Dispose of any used probe cover and store the unit in a safe place.
- NOTE:** Temperature readings vary from person to person and time to time; find your normal temperature and watch for changes. If you find large readings between two thermometers, it is possible to have different readings even both of them are within accuracy requirement.

RELIABLE METHOD TO MEASURE TEMPERATURE

THE TAKING OF BODY TEMPERATURES: Change in body temperature is one of the most important indicators of illness. It is important that accurate temperatures be taken. Even though the device is used in similar manner to the glass thermometers, it is important that the simple instruction for use of this device be followed closely and that the person taking the temperature understands body temperature.

ORAL USE: Place the probe well under the patient's tongue. Instruct the patient to keep their mouth closed while the thermometer is reading. A normal temperature by this method is usually considered between 36.0 °C to 37.5 °C.

RECTAL USE: If the rectal temperature is recommended by your doctor, you may use a probe cover. Insert the thermometer into the probe cover. Lubricate the probe cover with a water soluble jelly for easier insertion, do not use petroleum jelly, it may cause inaccuracy reading. Insert the tip of probe no more than 1/2 inch into rectum. STOP if you meet any resistance. The sensing unit is on the very tip of the probe and there is no need to insert the probe deep into the rectum. In general, the rectal temperature is around 0.3 °C to 0.5 °C higher than oral temperature.

NOTE: Due to petroleum is an oil soluble jelly, it will cause low temperature reading. If applying petroleum jelly is necessary, it is recommended to prolong the measurement regardless of the beep sound and continuous wait for around 2~3 minutes.

AUXILIARY USE: Wipe armpit with a dry towel. Place the probe in the patient's armpit and keep the patient's arm pressed firmly against the body. In general, the axillary temperature is around 0.5 °C to 1.0 °C lower than oral temperature. Note that in order to achieve better with axillary temperature measurement results, a minimum measuring time of 2~3 minutes is recommended regardless of the beep sound.

NOTE: Children should be attended by an adult during the entire measuring. After use, the digital thermometer should be stored in a safe place, out of reach of children. If you compare temperatures between two thermometers, it is possible to have different readings even if both thermometers are within their accuracy tolerance. Difference in how you take each reading and the thermometer's tolerance range can account for nearly 0.1 °C variability.

TROUBLESHOOTING

	Auto detection is erroneous. Please return to the distributor.
	Please redo measurements.
	The device detects low battery, please change the battery.

CLEAN & DISINFECTION

DO not use corrosive detergent to clean the unit. Clean and disinfect the thermometer before and after use.

- Cleaning process:
1. Clean the probe with cold water several times.
2. Clean the probe with a solution of 5~10% mild detergent about 10~15 minutes.
3. Continue to clean the probe with cold water until no chemical residues remain on the device.
4. Rinse the device with drinking water or distilled water to ensure the removal of all remnants of the cleaning agent.
5. If necessary, repeat step 2 through step 4 in order to make sure the device is clean.
6. A final wipe down by a clean wipe is necessary to avoid water stains.
- Disinfection process:**
1. To disinfect the device when completed the cleaning process.
2. Wiping the probe with a cloth dipped in a disinfectant solution such as 75% isopropyl (rubbing) alcohol : contact time $\geq$ 10 min for disinfection.
3. The disinfectant free time without residues.
- Drying process:**
1. After cleaning and disinfection process above.
2. Put the device to air-dry in the room temperature (20~30 °C), exposure time around 30~60 min.
3. And store the device in a dry, clean, and dust-free environment at room temperature.

CARE & STORAGE

- Do not drop or crush the device. It is not shock proof.
- Do not dismantle or modify the device other than changing the batteries.
- Store the device in an extreme environment (high or low temperature, high humidity or others) in a location with high humidity or severe dust. Please contact your local agents or distributors if the device does not function properly.

WARRANTY & SERVICE

The regular warranty period is 2 years and begins on the day of purchase. To benefit from an extended voluntary warranty period stated on the gift box, registration on our website is required. You can consult the full warranty terms as well as information on extending the warranty period and details of our services at www.bresser.de/warranty_terms.

LOT NUMBER INFORMATION

LOT NO. (MMYY) on the device label means the month and year of manufacturing. For example, "0324" means the device was manufactured in March, 2024.

BATTERY MANUFACTURER:

New Leader Battery
Industry (Deqing) Co., Ltd. Address: Chaoyang West Road No. 108, Deqing County, Zhaoqing City, P.R. CHINA
[www.newleader.com.hk](mailto:sales@newleader.com.hk) e-mail: sales@newleader.com.hk

BATTERY INFORMATION LR 41

Battery category: Button cell		
0% Mercury, Alkali-free, Zinc, Manganese.		
Chemical composition:	% Share in the composition	CAS Number
Manganese	30	1313-13-9
Graphite	3	7782-42-5
Potassium hydroxides	10	1313-58-3
Zinc	8	7440-66-6
Iron	37	7439-89-6
Water	11	7732-18-5
Paper	1	
Capacity/Weight: 43 mAh/0.63 g		
Fire extinguishing: Please extinguish with dry powder fire extinguisher. Wear a face mask.		

EC Declaration of Conformity

Bresser GmbH has issued a "Declaration of Conformity" in accordance with applicable guidelines and corresponding standards. The full text of the EU declaration of conformity is available at the following internet address: www.bresser.de/download/9810602/CE/9810602_CE.pdf

BATTERY REPLACEMENT

- When the symbol  appears, the battery needs to be replaced.
- Unhook the screw then pull the battery cover and the rubber plug off in the direction shown.
- Use a pointed object such as a toothpick or insulated probe to pry loose the old battery. Avoid using any sharp metal object in this operation.
- Place a new 1.5V button size battery, LR41 or equivalent in the chamber with the negative (-) side up and the positive (+) side down.
- Replace the rubber plug, the battery cover and the screw. Fasten the screw securely.

Warning:
1. Store battery and discard old battery carefully, out of reach of children. Swallowing the battery may be fatal. If the battery is swallowed, contact a hospital immediately to have it removed. Do not dispose of the battery in a fire. It may explode.

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POLYGREEN GERMANY GmbH